








## SURGICAL GOWN - WATERPROOF AND BREATHABLE - WASHABLE AT 95°C

Ref.: 570 800 02 00L000

### TECHNICAL SPECIFICATIONS

COMPOSITION		UNE EN ISO 1833:2011	
FABRIC	Polyester 100%	BASE:	Polyurethane 100%
STRUCTURE			
FABRIC	Knitted Interlock. Gauge 20. Polyester yarn 75/36/1 Deniers texturized . Total weight 65 g/m <sup>2</sup>		
BASE:	Waterproof & breathable membrane 15μ		
TOTAL WEIGHT	95 g/m <sup>2</sup> : 65 g polyester, 19 g polyurethane membrane and 11 g reactive polyurethane glue		
FINISH	Fabric with optical white and Hot Melt lamination with waterproof and breathable membrane		
CHARACTERISTICS			
DIMENSIONAL STABILITY	(Washed at 95°C) Weft -1%. Warp -2% UNE EN ISO 5077:2008		
WATERPROOFNESS	> 500 cm/H <sub>2</sub> O		UNE EN ISO 20811:1993
BREATHABILITY	> 150 g/m <sup>2</sup> .24h		BS 7209:1990
CERTIFICATIONS			
OEKOT-TEX® Class I ; OEKO-TEX® Class II ; 2018OK1127 AITEX / 2002AN4162 AITEX			
Company certified under ISO 9001:2015			
Velamen is member of SEDEX and audited under SMETA protocol ZC1058900			
AVAILABLE SIZES:		L	
AVAILABLE COLORS:		White	
CARE INSTRUCTIONS:		    	
30 WASHES AT 95°C GUARANTEED			
CONFECTION:			
Overlock elastic safety stitch. Six stitched straps to close the gown. Elastic rubber on the cuffs			



Backing: Mesh column structure



2 straps to adjust to the nape



4 ribbons to tie in pairs in inner and outer knots



elastic band at the ends of the sleeves to adjust on the wrists



Outer Side: polyurethane membrane laminated with the knitted polyester fabric

